
MIFU-014-02

Anaesthesia Reservoir Bags – Non Sterile
Medical device class; IIA

English

Devices used to provide a reservoir of gas for anaesthesia and ventilator machines. They connect to the machine using a 15mm or 22mm connection to BS EN ISO 5356-1.

Material: Polychloroprene – latex free

- 1) Do not use if packaging is damaged or wet
 - 2) To be used by personnel trained in anaesthesia/ventilation
 - 3) Ensure device is clear of any obstructions before use
 - 4) Ensure device is correct size for the patient
 - 5) Ensure the connector in the neck of the bag is secure inside the device
- Storage Temperature & humidity: 10°C (50°F) to 30°C (86°F), 30–50% RH
Operating Temperature: -18°C (0°F) to 50°C (122°F)

The devices are supplied non-sterile and are for single use or 7 days where appropriate filtration is added. Because of the risk of possible cross-contamination, these devices must not be reprocessed.

Discard the device after use in accordance with local regulations for medical waste. The device may be incinerated in a facility capable of burning plastic waste.

Adverse incidents may be reported to the Competent Authority in the country/region where the incident occurred.

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