

## ELECTROSURGERY POWER CABLES INSTRUCTIONS FOR USE

Picture 1 - Bipolar electro-surgical power cable (example)

Picture 2 Monopolar electro-surgical power cable (example)

Medical device Class 1S – sterile version      Class 1R – reusable version

### Intended use

1. Electro-surgical power cables are used to connect active electro-surgical devices to an energy generator. They cannot be dismantled.
2. Inspect the cable connection for damage or distortion.
3. Check the integrity of the machine end connection(s) by inserting them into the sockets of an electro-surgical generator or a suitably sized test socket. The plugs must provide a positive and secure connection to the machine. If in doubt withdraw the cable from service. When using a generator ensure that the generator is switched off before insert the cable plugs.
4. Check the integrity of the instrument connection by inserting a suitable instrument into the connection. The connector must provide a positive and secure connection to the instrument. If in doubt withdraw the cable from service. Before connecting the test instrument ensure that the cable is not connected to a power source.
5. Inspect the cable cover for cuts, abrasion or other damage. Look also for kinks in the cable especially at the points where the cable exits the moulded end connectors. Kinks and bulges in the cable cover, especially if the cable appears to be floppy in the area of the damage, indicates that the internal conductor may be damaged. If in doubt withdraw the cable from service.

### Intended user

Healthcare staff with appropriate training for electro-surgery procedures.

### Materials

- Multistrand copper conductor
- PVC insulation
- PP mouldings
- Brass / nickel plated plugs / sockets

THIS PRODUCT IS LATEX FREE.AND DOES NOT CONTAIN PHTHALATES

### Warnings

- This product should not be used if the integrity of the packaging has been compromised.
- Do not use if the packaging is contaminated or wet.
- This product is for **single use only**.
- The product must not be reprocessed.
- The product must not be altered or dismantled.

## ELECTROSURGERY POWER CABLES INSTRUCTIONS FOR USE

### Disposal

The devices may be disposed of using local procedures for contaminated medical waste.  
The materials are capable of being recycled if a safe process for recycling medical waste is available.

### Ordering information:

Meditech Systems Ltd, Unit 3 Shrublands Estate, Shearstock, Shaftesbury, Dorset, SP7 9PT, England

E mail; [orders@meditechsystems.co.uk](mailto:orders@meditechsystems.co.uk)

Tel: +44(0)1747 821546 Fax: +44(0)1747 825038

### EU Authorised Representative:

Clowrey Consultants, White House at Bridge, Holycross, Thurles, C Tipperary, Eire

E mail; [clowreyconsult@gmail.com](mailto:clowreyconsult@gmail.com) Tel; +44(0)7888744362

### ADVERSE INCIDENT REPORTING

Adverse incidents may be reported to the Competent Authority in the country / region where the incident occurred.

### Explanation of symbols;

	= Single Use. The definition of single use is for use on a single patient for a single procedure.		= Manufacturer
	= Do not reprocess / re-sterilise		= Temperature range
	= Use by date		= Consult instructions for use
	= Lot or batch number.		= Keep away from sunlight
	= Date of manufacture.		= Keep dry
	= Device catalogue number / product code / stock		= Contains latex
	= Caution		= Do not use if packaging is wet or damaged.
	= CE mark with notified body number		= Contains or Presence of Phthalates.

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