



## BERMAN INTUBATING / PHARYNGEAL AIRWAY

### Features:

- » Gold Standard difficult airway equipment
- » Large internal diameter accommodates the passage of an endotracheal tube through its center
- » May be opened on one side so it can be split and removed from around a tracheal tube, or similar.
- » Ideal for fibreoptic bronchoscopy and procedures for guiding, positioning and protecting expensive fibre optic equipment
- » Latex Free



**IT'S BACK!**

*Does your Difficult Airway Trolley need a Berman Intubating Airway?*

Product	Part No.	Quantity
Berman Intubating Airway, Size 8	22BIA008	ea
Berman Intubating Airway, Size 9	22BIA009	ea
Berman Intubating Airway, Size 10	22BIA010	ea
Berman Intubating Airway, Size 8 - Box of 10	22BIA008-10	Box of 10
Berman Intubating Airway, Size 9 - Box of 10	22BIA009-10	Box of 10
Berman Intubating Airway, Size 10 - Box of 10	22BIA010-10	Box of 10

Contact our team for further information and support

Meditech Systems Ltd

Unit 3 Shrublands Estate, Sherstock, Shaftesbury, Dorset, SP7 9PT  
Tel: +44(0)1747 821546 [www.meditechsystems.co.uk](http://www.meditechsystems.co.uk)