


CUSCO SPECULUMS



Strong construction with improved insulation performance for patient and user safety. Insulated and non-insulated options available. Smooth rounded design for improved patient comfort with a Secure centre self-retaining screw.





Main Features

- Single use high quality surgical steel instrument.
- Insulated & Non Insulated designs available.
- Nylon Insulated with Rilsan approved nylon.
- Secure centre screw design to ensure stable hold during use.
- Smooth blade design to ensure patient comfort.

Product	Size	Codes  (insulated)	Codes  (non-insulated)	Codes  (insulated)	Codes  (non-insulated)	Blade length(mm)	Blade length(mm)
---------	------	----------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------	------------------	------------------

GYNAECOLOGY INSTRUMENTS SPECULUM DESIGNS

Cusco	SMALL	32010*	32010N	SU32010	SU32010N	95	25
Cusco	MEDIUM	32011*	32011N	SU32011 	SU32011N	102	30
Cusco	LARGE	32012*	32012N	SU32012	SU32012N	110	35
Cusco	MEDIUM LONG	32013*	32013N	SU32013 	SU32013N	115	30
Cusco	EXTRA LARGE	32016*	32016N	SU32016	SU32016N	130	23

QUALITY | VALUE | INTEGRITY