

NASO PHARANGEAL



Naso-Pharangeal products are a single use, latex free product. We supply sterile which eliminates the risk of cross contamination. They are packaged in easy to open individual packs. Naso-pharangeals are made from PVC with DEHP and are available in blue colour. Available in 13 different sizes for comfortable patient fit.



Description	OD	ID	Code	Box qty
NASO PHARYNGEAL				
Naso pharangeal airway	14Fr	3.0mm	22507514B	10
Naso pharangeal airway	16Fr	3.5mm	22507516B	10
Naso pharangeal airway	18Fr	4.0mm	22507518B	10
Naso pharangeal airway	20Fr	4.5mm	22507520B	10
Naso pharangeal airway	22Fr	5.0mm	22507522B	10
Naso pharangeal airway	24Fr	5.5mm	22507524B	10
Naso pharangeal airway	26Fr	6.0mm	22507526B	10
Naso pharangeal airway	28Fr	6.5mm	22507528B	10
Naso pharangeal airway	30Fr	7.0mm	22507530B	10
Naso pharangeal airway	32Fr	7.5mm	22507532B	10
Naso pharangeal airway	34Fr	8.0mm	22507534B	10
Naso pharangeal airway	36Fr	8.5mm	22507536B	10
Naso pharangeal airway	38Fr	9.0mm	22507538B	10

QUALITY | VALUE | INTEGRITY