



BOUGIE/INTRODUCERS



Bougies allow positive location of the trachea and correct placement of the endotracheal tube (ET tube) in difficult airway situations. Bougies are a single use, latex free product. We supply sterile which eliminates the risk of cross contamination. They are packaged in easy to open individual packs. Bougies are made from Polyethylene (PE)



Clear Markings

The clear graduated marking provide depth indication reinforcing correct placement.

Reduced patient trauma

The rounded tip (coude) provides a smooth surface allowing the Bougie to be inserted without damaging the soft tissue, hence causing less trauma to the patient.

Coude Tip

The Coude tip allows for accurate placement by recognising the tracheal wall C-Rings. The flexibility allows for easy removal once the ET Tube is in place.

Description	FR	Length	Code	Box qty
-------------	----	--------	------	---------

BOUGIE/INTRODUCERS

Paediatric Bougie Coude Tip	10	70cm	220211PD	10
Paediatric Bougie Vented	10	70cm	220211PDV	10
Adult Bougie Coude Tip	15	70cm	220212AD	10
Adult Bougie Vented	15	70cm	220212ADV	10
Adult Bougie Vented	14	70cm	220212V14	10

QUALITY | VALUE | INTEGRITY